Thera-Band® Shoulder Pulley

Thera-Band products have been trusted by rehabilitation professionals for over 30 years. The enclosed shoulder pulley is specifically designed to assist patients undergoing shoulder rehab in regaining and maintaining range of motion. Patients can use the pulley during clinic visits or as part of a clinician-prescribed home exercise program.

Thera-Band Shoulder Pulley Features

- Distinct black marks on the rope for:
 - Immediate feedback on progress
 - Clear, accurate instructions
- A sturdy door anchor which is compatible with any standard door jamb
- Soft foam handles for superior comfort during use
- Adjustable rope length for custom fit and versatility

The Hygenic Corporation
Akron, Ohio 44310-2575 U.S.A.
(330) 633-8460
www.thera-band.com

For training examples go to www.Thera-BandAcademy.com

Thera-Band®, the Color Pyramid Design™, and the Associated Colors™ are trademarks owned by The Hygenic Corporation or its affiliates and may be registered in the United States and other countries. All rights reserved © 2010. Made in China. P04533

HOW TO USE THE MARKED ROPE

Reference to Specific Marks

The Thera-Band Shoulder Pulley features black marks along the rope. The black marks along either side represent 6 in/15 cm sections of rope. Each black mark can be referred to by number, with “#1” representing the first black mark from the center resting position, “#2” representing the second black mark from the center resting position, and so on. Please see the diagram to the right.

Immediate Feedback on Progress

At the start of each new exercise, make note of where the black marks fall relative to the pulley mechanism. Then make note of where the black marks are at the end point of the exercise. As range of motion improves, the black marks on the rope will provide immediate visual feedback on the progress made.

Clear and Accurate Instructions

In addition to providing visual feedback on progress, the black marks can be used for clinicians to provide clear, accurate instructions. Practitioners simply use the marks to set goals for patients and to provide guidelines so that patients do not push themselves too far too fast.


PRECAUTIONS

- The shoulder pulley is to be used for its intended purpose only.
- Before beginning any type of exercise program, consult your healthcare provider.
- **ONLY PERFORM EXERCISES PRESCRIBED BY YOUR THERAPIST OR PHYSICIAN.**
- This product is not a toy and should not be used without adult supervision. Consult your physician before allowing children to exercise.
- Before each exercise session, inspect the shoulder pulley for defects, cuts or wear. Replace the damaged product immediately.
- Consult your healthcare provider if you experience pain while performing your exercises.
- Hold the handles firmly to prevent them from slipping out of hands during use.
- Ensure that the connection to the door is secure prior to use.
- Do not wrap the rope of the shoulder pulley around any body part.
- Close and lock the door before starting your exercises. Make sure the door latches and locks and that no one opens the door while the pulley is in use.
- Ensure that the knot below each handle is secure prior to use.


GENERAL INSTRUCTIONS

- Use the Thera-Band Shoulder Pulley to perform exercises as instructed by a physician or therapist.
- Prior to beginning your exercises, you will need to secure your shoulder pulley to a door by following these simple steps:
 - Identify a door in your home that is not used often and that will allow you ample room to perform your exercises.
 - Place the Thera-Band webbing over the top of the door or along the side of the door, so that when the door is closed the round foam disc is on the opposite side of the door.
 - Close and lock the door before starting your exercises. Make sure the door latches and locks and that no one opens the door while the pulley is in use.
 - Avoid door gaps larger than 1” to maximize safety.
- The proper placement along the door jamb will be dependent on the goals of your prescribed exercises. Follow clinician guidance for any specific program.

Thera-Band Shoulder Pulley over the top of the door


Thera-Band Shoulder Pulley along the side of the door


- To adjust the length of the rope, tie a secure knot at the desired point below the handle.
- To monitor progress, track how close a given mark on the rope gets to the pulley mechanism during exercise.

Sample Exercise SHOULDER FLEXION


- Anchor the pulley to the top of the door
- Face the door and hold both pulley handles with a firm but relaxed grip
- Keep your hands shoulder width apart
- Pull down using your unaffected arm
- This will cause your affected arm to rise over your head
- Hold the stretch and repeat
- Use the black marks to monitor progress
- Face away from the door for extended range of motion

Sample Exercise SHOULDER ABDUCTION


- Anchor the pulley to the top of the door
- Face the door and hold both pulley handles with a firm but relaxed grip
- Keep the elbow of your unaffected arm close to your body and hold your affected arm out to the side
- Pull down using your unaffected arm
- This will cause your affected arm to rise in an arc
- Hold the stretch and repeat
- Use the black marks to monitor progress

Sample Exercise EXTERNAL ROTATION

- Anchor the pulley in the side of the door at elbow height
- Stand perpendicular to the door with your affected arm nearest the pulley
- Hold both pulley handles with a firm but relaxed grip, with your thumbs pointing up
- Place your arms against your sides and keep your elbows bent 90°
- Slowly pull outward with your unaffected arm, rotating your affected arm away from your body
- Make sure to keep the elbow of your affected arm firmly against your side at all times
- Hold the stretch and repeat
- Use the black marks to monitor progress

Sample Exercise INTERNAL ROTATION

- Anchor the pulley to the top of the door
- Hold the hand of your affected arm behind your hip (same side of body)
- Pull down gently using your unaffected arm
- This will cause your affected arm to move further up your back
- Hold the stretch and repeat
- Use the black marks to monitor progress

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